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STERILE FILLING NOZZLE for CBS™ high security straws

Intended use

The CBS™ filling nozzle is designed for the filling procedure of CBS™ High Security embryo and sperm straws for the prevention of contamination of the sealing area of the straw by the biological sample to guarantee airtight seals of the CBS™ High Security straw.

Cautions

- The user should read the entire Instructions for Use and understand the Cautions and Warnings and be trained in the correct procedure before using the CBS™ filling nozzle.

- Guarantee of the cryogenic resistance of the CBS™ High Security straws in liquid nitrogen depends on strict compliance with the instructions for use of the CBS™ filling nozzle and CBS™ High Security straws.

Warnings

- Manual filling of CBS™ High Security straws should be done with CBS™ filling nozzles in order to keep the sealing area clean and dry.

- When using CBS™ filling nozzles, never expel liquid back through the nozzle in order to avoid leakage between the filling nozzle and the straw and potential contamination of the sealing area.

- Only thermal seals made with SYMS type sealer with regular maintenance record (Cryo Bio System, France) are guaranteed air tight and cryopreservation resistant.

- Store in a clean and dry environment, away from light and heat sources.

Description

CBS™ filling nozzles are supplied as single use devices individually packaged in peel-off blister in a protective box per 100 units.

The CBS™ filling nozzle is composed of a transparent capillary overmolded with a flexible shape designed to connect airtight on the open end of CBS™ High Security straws. The capillary enters into the open end of the straw without coming into contact with the inside wall of the straw so as to keep a defined length of the straw clean and dry during the straw filling process and when withdrawing the filling nozzle. The small diameter and transparency of the capillary allows for selective aspiration of cells from a volume under microscope.

Filling and sealing


CBS™ High Security straws are filled by aspiration through a CBS™ filling nozzle. Aspiration should be done using a

micro-aspirator (Ref. 014498, Cryo Bio System.) or a syringe equipped with a special connector for CBS™ straws (Ref. 016730, Cryo Bio System).


1. Insert straight and firmly the side of the CBS™ High Security straw that is to contain the biological sample into the filling nozzle. A click can be perceived. The central capillary enters the straw without touching the inner surface. Check the connection by applying slight withdrawal force on the filling nozzle that should not move.

2. Load the biological sample into the straw by gentle aspiration force. Never expel liquid back through the nozzle in order to avoid leakage between the filling nozzle and the straw and potential contamination of the sealing area.

3. After completion of the filling of the straw, gently remove the filling nozzle using rotational

movements and avoiding any shaking and bending of the straw to avoid contact between the capillary and the inner surface. The two small openings on the side of the filling nozzle should not be blocked to enable air to enter when withdrawing the filling nozzle.

4. Visually inspect the sealing area for contamination before sealing the straw extremity with a SYMS type sealer. Disconnect the aspiration system and seal the second extremity of the straw as indicated above.

Disposal after use

Please follow specified legal local rules for disposal of biologically contaminated waste.

Quality assurance

Cryo Bio System is ISO 13485 certified.

CBS™ is a trademark of Cryo Bio System, IMV Technologies Group, FRANCE.

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