Ref. 022527

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Ref. 022527

CBS™ HIGH SECURITY TISSUE 2 ML STRAW

for cryopreservation
of tissue fragments

Caution

The user should read and understand the entire Instructions for Use and be trained in the correct procedure before using the CBS™ High Security tissue 2 mL straws.

Introduction

The physical properties and the method of filling and sealing CBS™ High Security straw ensure a leak proof environment after sealing for cryopreservation of biological samples of human origin, especially in liquid nitrogen.

Manufactured with biocompatible materials, the CBS™ High Security tissue 2 mL straw is designed for use in Assisted Reproduction Technologies, for example for cryopreservation of ovarian cortex.

Warning

Do not use if the package is opened or damaged as the sterility and the integrity of the product can no longer be guaranteed.

Store in a clean, dry environment, away from light and heat sources.

Cryo Bio System recommends to use the CBS™ High Security straws in a controlled environment, such as a laminar flow hood.

Description

The CBS™ High Security tissue 2 mL straw is a flexible, transparent tube made of non toxic ionomeric resin, 133 mm long and 5 mm in internal diameter. The CBS™ High Security straw should be thermally sealed with the equipment SYMS or SYMS II, purpose-designed by Cryo Bio System.

The CBS™ High Security tissue 2 mL straws are individually packaged and supplied sterile in bags of 20 units.

Identification

The CBS™ High Security tissue 2 mL straws can be identified by using LN2 resistant labels.

Manual filling and sealing

Caution: During filling the seal area of the straw should stay clean and dry.

1. Prepare the identification label and apply it on the outside of the CBS™ High Security tissue 2 mL straw, on the middle part of the straw at minimum 20 mm distance from the open extremes.

2. Holding the straw near the jaws, seal one extremity of the straw with the SYMS sealer equipped with the special guide for CBS™ High Security tissue 2 mL straws or directly with the SYMS II sealer.

3. To fill the straw with cryopreservation medium, use a disposable Pasteur pipette with diameter compatible with the diameter of the tissue straw. The operator must insert the pipette till the extremity of the CBS™ High Security straw and then

push out the medium slowly while pulling out the pipette at the same time. Watch out for any overflow, which are a potential source of contamination. It is important to use one new pipette each time there is a risk of contamination due to overflow.

4. Using a tong or very thin tweezers, carefully deposit the tissue fragment(s) into the straw, according to laboratory protocol. Tissue fragments should be small enough to enter the straw without spoiling the straw inside walls in the sealing area.

5. To assure dry and clean sealing area, wipe the inside of the extremity to be sealed with a sterile swab or pad.

6. Holding the straw near the jaws, seal the other extremity of the straw with the sealer SYMS or SYMS II.

Freezing

The freezing procedure of tissue fragments with the CBS™ High Security tissue 2 mL straw gives the same results with the added features of mechanical strength and tamperproof identification system (under the maximum conditions of security).

Storage in liquid nitrogen

A full line of storage containers provides versatility for the inventory management of the CBS™ High Security straws in liquid nitrogen: visotubes, goblets and canisters in different colors and/or dimensions are currently available.

Quality assurance

Cryo Bio System products are clean room manufactured from biocompatible materials and designed for medically assisted reproduction. Sterilization is done by irradiation.

Each lot of CBS™ High Security straws receives the following tests: endotoxin by LAL methodology and Mouse Embryo Assay.

All results are reported on a lot-specific certificate of analysis, which is available upon request.

Shelf-life

Three years from date of manufacture (see date mentioned on the label).

Disposal after use

Please follow specified legal local rules for disposal of contaminated medical waste.

CBS™ is a trademark of Cryo Bio System, IMV Technologies Group, FRANCE.

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