CS	Pokyny pro použití	2
	SLÁMKA CBS™ HIGH SECURITY NA SPERMA	
DA	Brugsanvisning	4
	CBS™ HØJSIKKERHEDSTRÅ TIL SÆD	
DE	Gebrauchsanweisung	6
	CBS™-HOCHSICHERHEITSPAILLETTE FÜR SPERMA	
EL	Οδηγίες χρήσης	9
	ΠΑΓΕΤΑ CBS™ ΥΨΗΛΗΣ ΑΣΦΑΛΕΙΑΣ ΓΙΑ ΣΠΕΡΜΑ	
EN	Instructions for use	11
	CBS™ HIGH SECURITY SPERM STRAW	
ES	Instrucciones de uso	13
	PAJUELA CBS™ DE ALTA SEGURIDAD ESPERMA	
ET	Instrucciones de uso	16
	CBS™ ÜLITURVALINE SPERMAKÕRS	
FI	Käyttöohjeet	18
	ERIKOISTURVALLINEN CBS™-SPERMAOLKI	
FR	Instructions d'utilisation	20
	PAILLETTE CBS™ HAUTE SÉCURITÉ SPERME	
HR	Upute za korištenje	23
	IZNIMNO SIGURNE SLAMČICE ZA SPERMU CBS™	
HU	Használati utasítás	25
	CBS™ HIGH SECURITY MŰSZALMA EMBERI SPERMA TÁROLÁSÁHOZ	
IT	Istruzioni per l'uso	27
	PAILLETTE CBS™ ALTA SICUREZZA SEME	
LT	Naudojimo instrukcija	30
	DIDELIO SAUGUMO CBS™ ŠIAUDELIS SPERMAI	
LV	Lietošanas instrukcija	32
	PAAUGSTINĀTAS DROŠĪBAS CBS™ SPERMAS SALMIŅI	
NL	Gebruiksaanwijzing	34
	PAILLETTE CBS™ HAUTE SÉCURITÉ SPERMA	
PL	Instrukcja obsługi	36
	SŁOMKA CBS™ O WYSOKIM POZIOMIE BEZPIECZEŃSTWA DO NASIENIA	
PT	Instruções de utilização	39
	Palheta CBS™ de Alta Segurança sêmen	
SK	Návod na použitie	41
	VYSOKO BEZPEČNÁ SPERMIOVÁ BANKA CBS™	
SL	Navodila za uporabo	44
	SLAMICA CBS™ Z VISOKO STOPNJO VARNOSTI ZA SEMENSKO TEKOČINO	

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EN

CBS™ HIGH SECURITY SPERM STRAW

Cryopreservation device for sperm

Indications

CBS™ High Security sperm straws are used for cryogenic slow freezing and liquid nitrogen storage of human sperm samples in assisted reproduction procedures.

Intended use

The CBS™ High Security sperm straw is specifically designed for the cryopreservation of human biological samples in liquid nitrogen including sperm.

Clinical benefit

Cryo bio system does not claim any clinical benefit for CBS™ high security straws, only the indirect benefit of enabling an infertile couple to conceive

Contraindications and side effects

There are no known contraindications or side effects to the use of this device.

Target population

Patients suffering from a dysfunction of their reproductive system that prevents them from conceiving by natural means.

Patients freezing their gametes prior to treatment of a serious pathology that could affect their ability to conceive through natural means.

Persons wishing to cryopreserve sperm for the purpose of medically-assisted reproduction.

Users

The personnel using this device must be qualified, especially regarding ART laboratory Best Practices and nitrogen handling according to regulations in force.

Warnings

A sterile filling nozzle for CBS™ High Security straw (Ref. 024856, Croy Bio System) must be attached to the CBS™ High Security straw before (manual) filling. Be sure to follow the instructions for use of the nozzle for its connection and disconnection.

Never push back media after applying suction so as to prevent the passage between the filling plug and the straw, thus preventing the humidification of the sealing area.

A characteristic tactile and/or auditory «click» crimps the straw correctly to the filling nozzle.

The CBS™ High Security sperm straw is a single-use sterile device.

This product cannot be reused, as any reuse of a product is likely to entail risks of cross-contamination between different patients.

Do not use the product if the packaging is damaged or open as the sterility and integrity of the product can no longer be guaranteed.

Do not use a product that has passed the expiry date on the packaging.

Do not resterilize the product, even if the packaging is damaged. Re-sterilization may contaminate the straw and the biological sample and affect the mechanical properties of the straw.

Only seals made with a SYMS sealing machine (Cryo Bio System, France) are guaranteed to be leakproof and resistant to cryopreservation.

If the device malfunctions, it must no longer be used and must be conserved as part of a potential investigation. A new device must be used.

Precautionary measures

The user should read the entire Instructions for Use before using the straw.

Store in a clean, dry place, away from heat sources.

Description

The physical properties and method of sealing, CBS™ high security sperm straws guarantee that they remain sealed when immersed and stored in liquid nitrogen, with no direct contact between the biological sample and liquid nitrogen.

CBS™ high security sperm straws are for single usage and are packed in 4 or 5 individual peel-off blisters providing a sterile barrier, which are placed in an outer bag and sterilized by irradiation.

The CBS™ High Security sperm straw has a flexible, clear, non-toxic ionomer resin tube, with length 133 mm and inner diameter 2.5 mm, which contains a braided safety plug. It can be autogenically sealed using the SYMS type sealer, PACE, MAPI or DIVA equipment, specifically designed for this purpose by Cryo Bio System:

CBS™ 0.3 mL model

This straw comprises two separate parts separated by a white plug.

The first section is designed to hold the sperm sample and has a useful volume of 0.3 mL.

The second section, on the opposite side of the useful volume, allows for manipulation and identification of the straw.

For use with PACE or MAPI, CBS™ High Security 0.3 mL straws are provided with fixed colored identification rods and white external jacket.

CBS™ 0.5 mL model

The section designed to hold the sperm sample has a useful volume of 0.5 mL. The external part of the sliding security plug is colored.

Equipment required but not included:

Sealing machine (only heat sealing machines from the SYMS® range have been approved as compatible).

Micropipetter (ref. 014498, Cryo Bio System) or a syringe fitted with a specific connection tip (ref. 016731 or 016730).

Filling nozzle (ref. 024856, Cryo Bio System).

Automatic filling, sealing and printing machine (MAPI, ref. 007218).

Semi-automatic filling and sealing machine (PACE, ref. 018451).

Injection nozzles (ref. 007454 or ref. 018620).

aspiration nozzles (ref. 007451).

Quality

Ability to be filled correctly

Ability to be printed and sealed

Ability to withstand nitrogen

Each batch of CBS™ high security sperm straws undergoes the following tests:

Determination of endotoxins using LAL method (<20 EU/unit)

MEA test (> 80% expansion rate to blastocysts in 96 hours).

A- Identification of CBS™ High Security Sperm Straws

The biological sample must be identified according to laboratory protocols.

Depending on the straw reference used, it is possible to:

Identify by printing (with MAPI machine) or manually identify (with a nitrogen-resistant ink pen) on the white jacket.

Identify by sticking a nitrogen-resistant label on an identification rod which will be introduced in the straw on the unused side after filling. After sealing the straw, the identification is then tamper-proof. (step B.5).

Identify by sticking a nitrogen-resistant label on the outside of the straw (without jacket).

B- Manual Filling and Sealing of CBS™ High Security Sperm Straws

For sperm storage with a view to specific applications such as intrauterine insemination or in-vitro fertilization (IVF), CBS™ High Security sperm straws are filled individually using a specific sterile filling nozzle and a micropipetter or a syringe equipped with a special connector. The SYMS sealer assures an impermeable seal.

- 1- Attach a sterile filling nozzle to the straw at the opposite end from the plug. A characteristic click can be felt and/or heard when the nozzle is correctly fitted.

- 2- Attach a micropipetter or a syringe equipped with a special connector on the plug end. Aspirate the prepared sperm. The biological sample should rise through the first cotton and up to half of the powder.

3- Gently withdraw the filling nozzle by twisting it, avoiding any shaking movements. This is to keep the end of the straw free of sperm.

4- Seal the released extremity with a SYMS type sealer, holding the straw near the sealing jaws.

5- Seal the other extremity as indicated above with a SYMS type sealer.

B'- Automatic filling and sealing of CBS™ High Security Sperm Straws

For large-scale application (for example, sperm banks), the CBS™ High Security sperm straws are filled with PACE, MAPI or DIVA automated filling and sealing machines using a purpose-designed injection and aspiration nozzle. Refer to the user manual of the equipment used.

C- Freezing CBS™ high security straws

Once the straw is filled and sealed on both sides, freeze according to laboratory regulations and refer to the user manual of the equipment used.

Personnel must be qualified in nitrogen handling as per the applicable legislation in the country.

D- Storage in liquid nitrogen

CBS™ High Security sperm straws are protected from the mechanical constraints of being stored in liquid nitrogen by using CBS™ visotubes, goblets and canisters.

The cold chain must always be maintained when carrying out straw freezing, storage and identification procedures in order to ensure the integrity of the biological sample

E- Thawing

To prevent any risk of unwanted thawing, ensure that the part of the straw containing the biological sample remains in the liquid nitrogen at all times

1. Verify the identification of the straw to be thawed (by the color of the internal identification rod and/or its alphanumeric and/or bar code).
2. Carefully extract the selected straw from the visotube in which it is stored with a small forceps and transfer it immediately to a water bath (+37°C).
3. Thaw the straw by immersion in warm water (+37°C maximum) for about 30 seconds.
4. Make sure to proceed with the thawing one straw at a time in order to preserve traceability.

Usage in Intra-uterine insemination / In-vitro fertilization

1. Prepare the working area and the receptacle (Tube, Petri Dish) previously identified according to the straw, based on the destination of the biological sample and the protocol used in the laboratory.
2. After having wiped and decontaminated the straw, empty it by capillarity action on top of the receptacle. For this:
Cut with a dry single-use bistouri or decontaminated scissors the distal sealed part with respect to the plug, then partially cut or open the straw just below the plug to avoid cutting inside the plug.

Disposal after use

Please follow specified local regulations for disposal of contaminated medical waste.

Reporting of serious incidents

Any incident related to the device should be reported to CBS at the following address materiovigilance@imv-technologies.com. Serious incidents must be reported to the competent authority of the Member State in which the user is established.

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