

# High Security Vitrification

## Integral HSV

Improved cooling speed  
Uninterrupted cold chain



- Vitrification
- Vertical sealing
- Storage area
- Waiting area
- Opening area

**Connected to**...LIFE SCIENCES


**Cryo Bio System**  
Groupe I.M.V. Technologies

# High Security Vitrification


## Integral HSV




Removable and articulated daisy goblet basket ensuring an unbreakable cold chain



Waiting area for straws before opening and warming-up



LN2 level indicator



Storage area for accessories

Straw opening area



Vitrification slots connected to color-coded and patient-identified straw holders



**Connected to... LIFE SCIENCES**

[www.cryobiosystem-imv.com](http://www.cryobiosystem-imv.com)