



CBS™ HIGH SECURITY VITRIFICATION STRAWS

Exclusive High Security closed system for embryo vitrification



For the highest level of safety and quality in vitrification


A circular graphic composed of ten colored segments in a ring, with colors including blue, purple, pink, red, orange, yellow, green, and light blue. The segments are arranged in a circle around a white center.

Medical Device Class IIa

CBS™ HIGH SECURITY VITRIFICATION STRAW

Cryo Bio System

1. High Security ultra-thin outer straw
Pre-sealed weighted end and flared opening
2. Capillary tube with pre-formed fully clear gutter
co-injected with its colored identification rod
3. Blue introduction device



SYMS III



Improved performance with the use of our **Integral HSV**



• RELIABILITY

Thermo-sealed ionomeric resin straw

• TRACEABILITY

Color coding and LN2 resistant labels for easy storage and identification

• SECURITY

Straightforward simple use

• EFFICIENCY

Safe and fast procedure

• EASY STORAGE

Compatible with CBS™ storage concept

• ONE-STEP OPENING AND RETRIEVAL



INTERNATIONAL STANDARDS COMPLIANCE

CE, FDA 510(k) cleared



+33 (0) 233 346 464



www.cryobiosystem.com



contact@cryobiosystem-immv.com



ZI n°1 Est, 61 300 L'Aigle, FRANCE